

SCIENCE IMAGING SYSTEMS

Application Note

 No. 7

Film Method and CCD Method for Chemiluminescence Detection

LAS-1000

Introduction

The chemiluminescence method of detecting DNA or protein on blotted membrane filters is widely used. Until recently, this method depended almost exclusively on x-ray film for detection of the chemiluminescent emission. This is now being changed by dramatic advances in cooled CCD camera technology. Today, the traditional x-ray film system is gradually giving way to high-sensitivity cooled CCD camera systems.

Both systems have the same purpose - to detect the emitted light - but they differ totally in principle. The film method records the pattern of the light by a chemical reaction while the CCD method converts the light into an electrical signal. The difference in principle leads to different image features. The CCD method is superior in a number of points: the digital image it produces is broader in dynamic range, has better linearity and therefore more accurate.

This Application Note compares the light detection principles of the film method and the CCD method and explains the imaging processes.

Contents

1. Light Detection and Imaging
2. Comparison of the Film Method and CCD Method

Summary

- The film and CCD methods are the ones generally used for chemiluminescent imaging.
- The digital image obtained by the CCD method is ideal for quantitative analysis.
- With the CCD method, imaging is simple and optimum exposure conditions are easy to determine.
- When equipped with the URF20L (f0.85) lens, the LAS-1000 provides sensitivity equal to, or better than, the film method.

1 Light Detection and Imaging

Film Method

The chemiluminescent light forms an image on x-ray film by reducing the silver halide grains in the film's emulsion layers followed by development procedure.

(1) Layer structure of the film




Fig.1-1

(2) Film imaging process




Fig.1-2

Exposure to the image-carrying light generates development centers in the emulsion layers in the pattern of the image. Development produces a visible image by reduction of the silver halides near the development centers to silver metal.

CCD Method

The CCD camera has a matrix of over one million light-receiving elements each corresponding to one image pixel. The elements (photo pixels) accumulate electrons in proportion to received light and then convert the electrons to an electrical signal.

(1) CCD photo element




Fig.1-3

When light enters a CCD element, the element generates and accumulates electrons.




Fig.1-4

Silver halide

A compound of silver (Ag) and a halogen. AgCl, AgBr and AgI are typical silver halides. General symbol: AgX.

Emulsion

Short for "photographic emulsion," a suspension of photosensitive silver halide grains in an aqueous gelatin solution.

Fig.1-1 Simplified diagram of x-ray film layers

Development

To obtain a visible image. Silver halides are reduced to silver metal by a reducing compound such as phenidone or hydroquinone.

Fixing

For stabilizing the developed image. A weak acidic solution of sodium thiosulfate is used to dissolve undeveloped silver halide grains.

Fig.1-2 Schematic diagram of the film imaging process

Gate

Voltage is applied to this electrode to transfer the accumulated electrons as image data.

Oxide film

A thin, transparent layer of silicon dioxide for insulating the gate from the substrate.

Silicon substrate

Most CCD image sensors are fabricated on a base of silicon crystal.

Fig.1-3 Simplified, conceptual diagram of CCD element (photo pixel)

Fig.1-4 CCD imaging process

The electrons accumulated on the CCD elements (photo pixels) are sent to an output section through transfer lines and then used as image data by a computer.

2 Comparison of the Film Method and CCD Method

Light Emission and Image-capture

In the film method, the sample is placed in direct contact with the film. This maximizes the exposure to light, however not all of the photons entering the film are used to reduce the halide grains of silver. For instance, several photons are needed to form a development center.

In the CCD method, light from the sample is collected through a lens. This decreases the light emission entering the CCD. However, the CCD is very efficient. That it generates and stores an electron for almost every photon received.

Obviously, the detection sensitivity of the CCD method depends greatly on how effectively the light is collected. That is why the high-sensitivity Fujinon URF20L (f0.85) lens has been made available for use with the LAS-1000.

Sensitivity

The LAS-1000 was used to detect extremely weak chemiluminescence [photon flux: approx. 5×10^5 photons/sec \cdot cm 2 (equivalent to 5 fg spot in the below image.)]. These results are comparable to those obtained with film in the same concept.




Fig.2-1

Imaging Process

(1) Film method




Fig.2-2

In the film method, the steps from exposure through development are conducted in a darkroom.

(2) CCD method





Fig.2-3

In the CCD method, a digital image can be obtained in a light room by conducting the exposure in a special, light-tight darkbox like that in the LAS-1000.

Film method exposure



CCD method exposure




Fig.2-1 Detection of chemiluminescence by film method and CCD method

Sample : pBR328 dot blotted membrane filter
 Reagent : CDP-Star®
 Film : FUJI newRX
 Detection system : LAS-1000
 Lens : URF20L (f0.85)

Buffer was applied as negative control. (0fg)

Fig.2-2 X-ray film imaging process

Fig.2-3 CCD imaging process

■ Quantitative Image Analysis




Fig.2-4

Fig.2-4 Light intensity vs image intensity with x-ray film and CCD

The density of the image by the film method is not directly proportional to the light intensity and only a segment of the curve can be treated as linear.

The intensity of the image by the CCD method is directly proportional to the light intensity. Quantitative analysis is greatly simplified.

■ Sample Size

The largest sample that can be imaged by the LAS-1000 is determined by the tray height. Since the CCD has fixed 1384 x 922 available pixels, the resolution of the image is determined by the tray height. LAS-1000 is designed to get enough resolution for practical recognition of separate bands on the blotted membrane.

Tray height vs. max. sample size in LAS-1000

Tray height	Max. sample size (cm)
1	7 x 11
2	9 x 14
3	12 x 18
4	14 x 21
5	17 x 25
6	20 x 25
7	25 x 25

Tray height

The sample tray of the LAS-1000 can be set at any of seven levels.

Editors & Writers

Kenji Miura, Ph.D.

Tohru Tsuchiya, M.Eng.

Makiko Nagashima

Chiori Imai

(Fuji Photo Film Co., Ltd.)

CDP-Star is a registered trademark of Tropix Inc.

March 1998



FUJI PHOTO FILM CO., LTD.

SCIENCE SYSTEMS, EQUIPMENT PRODUCT DIVISION

26-30, NISIAZABU 2-CHOME, MINATO-KU, TOKYO 106-8620, JAPAN

Telephone: +81-3-3406-2201

Facsimile: +81-3-3406-2158

E-mail: sginfo@tokyo.fujifilm.co.jp

©1997-1998 Fuji Photo Film Co., Ltd.